**Tedequim**

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**TEDEQUIM® OHCal Hidróxido de Calcio Polvo**

**PÁGINA..... 2; 3**

**INSTRUCTIONS FOR USE**

**TEDEQUIM® OHCal**

**CALCIUM HYDROXIDE POWDER - PA**

**PAGE..... 4; 5**

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**Rev. B**

## **OHCaI CALCIUM HYDROXIDE POWDER - PA**

Calcium hydroxide is considered the best medication to induce the deposit of hard tissue and favor the healing of the vital pulp and periapical tissues.

Both the biological and clinical properties of calcium hydroxide,  $\text{Ca(OH)}_2$ , are associated to its physicochemical properties. The first are related to the ionic activity, due to the dissociation in the ions  $\text{Ca}^{++}$  and  $\text{OH}^-$  and their high pH since it is a base, and this confers its bactericidal property. The second, to the low solubility of this base. The water solubility of  $\text{Ca(OH)}_2$  is 0.8 g / L at 25 °C and decreases with increasing temperature. This low solubility is a good clinical characteristic as it takes a long time before it is solubilized in the tissue fluid when it is in direct contact with the vital tissue.

This rate of solubility or reabsorption will also depend on the vehicle in which the calcium hydroxide is set.

For endodontics use it is convenient that the paste is formed by the calcium hydroxide powder, a radiopaque material and an appropriate vehicle.

Calcium hydroxide pastes are classified into three categories depending on the vehicle used. The aqueous, which have water as a vehicle; viscose, whose vehicle are water-soluble substances; and the only ones that have water insoluble compounds such as oils and fats. Each of them have different clinical properties.

In summary, in those clinical situations that require a rapid ionic release at the beginning of the treatment, a paste of calcium hydroxide with aqueous vehicle will be more suitable while, in clinical situations that require a gradual and uniform release of ions, the viscous vehicle is more convenient. Finally, pastes with oil vehicle have more restricted use since they will only be used when a very slow ionic dissociation is required.

### **Dental use of Calcium hydroxide is:**

- As a coating in large cavities with infected dentin when the state of the pulp is uncertain and the dentin is sensitive due to the absence of irregular secondary dentin. Also in hypoplastic areas or fractured incisors.
- As an antiseptic.
- For direct and indirect pulp coating.
- To try the induction of a continuous apical closure.
- Temporary filling after pulpectomy, protective filling of the root canals after the instrumentation, washing and proper treatment of the walls of the root canals.
- To induce the formation of repair dentin.
- The powder must be mixed with distilled water to form a paste, whose pH is 12.5, which would fulfill the function of the aqueous OHCaI paste.

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- Temporary filling of root canals.
- Intermediate fillings between endodontic sessions.
- In all clinical situations where alkalization and rapid diffusion of calcium hydroxide should be required.
- Partial biopulpectomy

**Instructions for use:**

1. Remove the calcium hydroxide powder using the dispenser and place it in a suitable container. Add distilled water (or glycerol) and mix until a paste with the required consistency is obtained.
2. Use TEDEQUIM EDTA solution to remove.
3. Regardless of the indication and the presentation used, it must be remembered that the properties of Ca hydroxide are effective only when it comes into direct contact with the area to be treated. Therefore, care should be taken to apply and condensate the applied material; that can be done with apical condensers, probes soaked in cotton, or big caliber paper cones.

**Storage:**

Store in a cool, dry place. Keep between 4° and 30° C.

**Authorized by ANMAT PM 1277-23**

**Technical Director: Miryam Della Vedova - Pharmacist - License 2107**

**Exclusive sale to professionals and sanitary facilities**

**Manufactured by TEDEQUIM SRL**

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**X5147GGP, Córdoba, Argentina.**

**Tel: +54-3543-448260**

**[www.tedequim.com.ar](http://www.tedequim.com.ar)**

**Rev. B**